

Non-invasive Measurement of Intestinal Fat Absorption C1059

Version: 1

Edited by: Patrick Tso, Dana Lee

Summary
Reagents and Materials
Protocol

Summary:

Dietary fat containing 5% sucrose polybehenate, a non-absorbable food additive, is fed in a semi-synthetic diet containing absorbable fat to rodents. The animals eat the diet ad lib for 3 days. Two or more fecal pellets are collected on days 3 and 4. Fat absorption is calculated from the ratios of behenic acid to other fatty acids in the diet and in the feces as analyzed by gas chromatography of fatty acid methyl esters.

Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Rodent Diet With 45 kcal% Fat and Modification To 16% Fat	Research Diets	D07031501-02
(Mostly Safflower Oil or Butter Fat) With 5% Fat as Sucrose		
Polybehenate		

Protocol:

Animal Feeding and Fecal Collection:

MONDAY: REPLACE CHOW WITH TEST DIET

TUESDAY: CONTINUE TEST DIET,

NEW CAGE BEDDING BEFORE FEEDING PERIOD

WEDNESDAY: CONTINUE TEST DIET (DAY 3 – RETURN)

COLLECT 5-8 FECAL PELLETS/ANIMAL

CHANGE CAGE BEDDING BEFORE FEEDING PERIOD

THURSDAY: COLLECT 5-8 FECAL PELLETS/ANIMAL

(DAY 4 – RETURN)