



# D3406 HDL-C and LDL-C/VLDL-C

And

# D3407 to "HDL-TG and LDL-TG/VLDL-TG Protocols

Version: 1

Replaced by version N/A

Edited by: Peter Havel - UC Davis Metabolism and Endocrinology Core

[Summary](#)

[Reagents and Materials](#)

[Protocol](#)

[Reagent Preparation](#)

## Summary:

LDL and VLDL are separated from HDL using a precipitation reagent. Then the HDL fraction is measured for either TC or TG using the same reagents for total cholesterol or triglyceride.

## Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Calibrator	Fisher Diagnostics	TR43002
TC Reagents	Fisher Diagnostics	TR13421
TG Reagents	Fisher Diagnostics	TR22421
2X LDL/VLDL Precipitation Buffer	AbCam	ab105138
PBS		
Microplate		
Platereader		

## Protocol:

1. Add 25µl 2X precipitation buffer to 25µl of sample using a positive displacement pipet.
2. Vortex and let sit at RT for 10 minutes.

3. Centrifuge at 2000×g for 10 minutes at 4°C.
4. Pipet supernatant into new tube, this is the HDL fraction.
5. Add 5 µl of calibrator and sample to each well.

*IMPORTANT: Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.*

6. Add 300 µl of TC or TG reagent to each well. Incubate at 37°C for 5 minutes. Read at 540 nm.

*IMPORTANT: If samples are hemolyzed, pipet a blank well with 5µl sample and 300µl PBS*

7. Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as  $(\text{Sample Absorbance} \div \text{Calibrator Absorbance}) \times \text{Calibrator Concentration}$ .
8. HDL samples are diluted 1/2 so multiply these by 2 to get the final value. Subtract this from the total triglyceride or cholesterol value to get the LDL/VLDL value.

## Reagent Preparation:

Reagent – ready to use

PBS – ready to use