



# NEFA-HR Assay

## C1057

Version: 1

Edited by: Patrick Tso, Dana Lee

[Summary](#)

[Reagents and Materials](#)

[Protocol](#)

[Reagent Preparation](#)

[Working Color Reagent Solutions A](#)

[Working Color Reagent Solutions B](#)

### Summary:

Quantitative determinations of non-esterified fatty acids in plasma/serum/lymph will be made using the NEFA-HR enzymatic colorimetric method assay.

### Reagents and Materials:

Reagent/Material	Vendor	Stock Number
HR Series NEFA-HR(2) Color Reagent A	Wako	999-34691
HR Series NEFA-HR(2) Solvent A	Wako	995-34791
HR Series NEFA-HR(2) Color Reagent B	Wako	991-34891
HR Series NEFA-HR(2) Solvent B	Wako	993-35191
NEFA Standard Solution	Wako	276-76491

### Protocol:

- Prepare working Color Reagent Solutions A and B.
  - Reconstitute **Color Reagent A** with a portion of **Solvent A** and then transfer entire contents into **Solvent A** bottle, rinsing Color Reagent vial several times.
  - Reconstitute **Color Reagent B** with a portion of **Solvent B** and then transfer entire contents into **Solvent B** bottle, rinsing Color Reagent vial several times.
- Locate working Standard (1mmol/L or 1 mEq/L).  
**THIS ASSAY DOES NOT REQUIRE A SERIAL DILUTION**
- Using a 96 well flat bottom plate, into separate wells, pipette 5 $\mu$ L of deionized water, 1mM standard, or sample to be assayed.
- Add 200 $\mu$ L of **Color Reagent Solution A** to all wells.
- Mix well and Incubate plate for 5 minutes at 37°C.
- Measure the absorbance of each well at 550nm (sub:660nm). This measurement (Abs1) will serve as the sample blank.
- Add 100 $\mu$ L of **Color Reagent Solution B** to all wells.
- Mix well and Incubate plate for 5 minutes at 37°C.
- Measure the absorbance of each well at 550nm (sub:660nm). This will be your Abs2 value.
- Obtain the final absorbance (Sample<sub>abs</sub>) by subtracting the first reading (step 5) from the second reading (step 8).\*

10/17/13

11. Plot the absorbance vs. concentration to construct the calibration curve. A linear calculation model should be used.
12. To calculate sample concentration by calculation use the following formula:

$$\text{Sample Conc.} = (\text{Sample Absorbance}/\text{Standard Absorbance}) * \text{Standard Concentration}$$

\*The sample blank absorbance (Abs1) from the first measurement (step 5) should be multiplied by a Factor (F) in order to correct for changes in volume, as follows:

$$F = (\text{Sample vol} + \text{R1 vol}) / (\text{Sample vol.} + \text{R1 vol} + \text{R2 vol})$$

$$\text{For this assay: } F = (5+200) / (5+200+100) = \mathbf{0.67}$$

$$\text{Therefore: } \mathbf{\text{Sample}_{abs} = \text{Abs2} - (\text{Abs1} * 0.67)}$$

**Specimen:** Serum or Plasma. Specimen stable for 7 days at 2-8°C or 3 months at -20°C.

**Assay Linearity:** 4.0 mEq/L

**Reagent Stability:** 7 days at 2-8°C

**Stability of Final Reaction:** 60 minutes

## Reagent Preparation:

[Working Color Reagent Solutions A](#)

[Working Color Reagent Solutions B](#)

### Working Color Reagent Solutions A:

Reagents and Materials:

**Color Reagent A**

**Solvent A**

Procedure:

Reconstitute **Color Reagent A** with a portion of **Solvent A** and then transfer entire contents into **Solvent A** bottle, rinsing **Color Reagent** vial several times.

### Working Color Reagent Solutions B:

Reagents and Materials:

**Color Reagent B**

**Solvent B**

Procedure:

Reconstitute **Color Reagent B** with a portion of **Solvent B** and then transfer entire contents into **Solvent B** bottle, rinsing **Color Reagent** vial several times.