



# Phospholipids Assay

## C1060

Version: 1  
 Edited by: Patrick Tso, Dana Lee

- [Summary](#)
- [Reagents and Materials](#)
- [Protocol](#)
- [Reagent Preparation](#)
- [Working Reagent](#)

### Summary:


Determinations of phospholipids in plasma/serum/lymph will be made using the Wako Phospholipids C enzymatic assay. In this assay, phospholipids in the sample are hydrolyzed ultimately producing a blue pigment. The amount of phospholipids in the sample is determined by measuring the absorbance of the blue color.

### Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Phospholipid C assay kits, standard included in kit	Wako	433-36201

### Protocol:

1. Prepare working standards by making a serial dilution of the stock 300mg/dl standard. (See diagram below) \*Stock standard included in kit.



2. Prepare **Working Reagent** by reconstituting one vial of **Color Reagent** with a portion of **Buffer** then transferring entire contents to **Buffer** bottle, rinsing **Color Reagent** vial several times.

3. Using a 96 well flat bottom plate, into separate wells, pipette 2 $\mu$ L of deionized water, standard, or sample to be assayed.
4. Add 300 $\mu$ L of **Working Reagent** to all wells.
5. Incubate plate for 5 minutes at 37°C.
6. Determine the absorbance (abs) of the standards and of each unknown at 600nm.
7. Calculate values of unknowns from the standard curve.

**Specimen:** Serum or Plasma. Specimen stable for 7 days at 2-8°C or 3 months at -20°C.

**Assay Linearity:** 1000 mg/dl

**Reagent Stability:** 7 days at 2-8°C

**Stability of Final Reaction:** 60 minutes

## Reagent Preparation:

### [Working Reagent](#)

#### **Working Reagent:**

##### Reagents and Materials:

**Color Reagent**

**Buffer**

##### Procedure:

Reconstitute one vial of **Color Reagent** with a portion of **Buffer** then transferring entire contents to **Buffer** bottle, rinsing **Color** vial several times.