

# Triglyceride Assay C1092

Version: 1  
 Edited by: Patrick Tso, Dana Lee

- [Summary](#)
- [Reagents and Materials](#)
- [Protocol](#)
- [Reagent Preparation](#)
- [Working Reagent](#)

## Summary:


Determinations of triglycerides in plasma/serum/lymph will be made using a Randox Triglycerides colorimetric kit. The triglycerides are determined after enzymatic hydrolysis with lipases.

## Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Triglyceride Assay Kit	Randox	TR213

## Protocol:

1. Prepare **Working Standards** by making a serial dilution of the stock 200mg/dl standard. (See diagram below) \*Stock standard is included in kit



2. Prepare **Working Reagent** by reconstituting one vial of **Enzyme Reagent R1b** with a portion of **Buffer R1a** and then transfer entire contents to bottle **R1a**, rinsing vial **R1b** several times.
3. Using a 96 well flat bottom plate, into separate wells, pipette 2µL of deionized water, standard, or sample to be assayed.

4. Add 200 $\mu$ L of reconstituted **Working Reagent** to all wells.
5. Incubate plate for 5 minutes at 37°C
6. Determine the absorbance (abs) of the standards and of each unknown at 500nm.
7. Calculate values of unknowns from the standard curve.

**Specimen:** Serum or Plasma. Specimen stable for 7 days at 2-8°C or 3 months at -20°C.

**Assay Linearity:** 1172 mg/dl

**Working Reagent Stability After Reconstitution:** 21 days stored at 2-8°C. PROTECT FROM LIGHT

**Stability of Final Reaction:** 60 minutes

## Reagent Preparation:

[Working Reagent](#)

### Working Reagent:

Reagents and Materials:

**Enzyme Reagent R1b**

**Buffer R1a**

Procedure:

Prepare **Working Reagent** by reconstituting one vial of **Enzyme Reagent R1b** with a portion of **Buffer R1a** and then transfer entire contents to bottle **R1a**, rinsing vial **R1b** several times.