



Creatinine

Version: 1

Edited by: Gary Cline, John Stack: Yale MMPC Analytical Core

Summary: Procedure followed to detect the concentration of creatinine in serum, plasma, and urine.

Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Acetonitrile	J.T. Baker	9012-03
LC/MS/MS Buffer		
D3-Creatinine	CDN Isotopes	W212P16

Protocol:

1. Fill 1.5mL eppendorf tubes with 80 μ L of acetonitrile.
2. Add sample:
 - a. Use 20 μ L of plasma or serum/tube
 - b. Use 6 μ L of urine/tube
3. Add 6.5 μ L of d3-creatinine (5mg/dl) to each tube.
4. Vortex each sample for 5 seconds.
5. Spin each sample at 4°C for 10 minutes
6. Transfer supernatant to LC/MS/MS vials.
7. Dry samples in speed vacuum.
8. Resuspend pellet in 75 μ L of LC/MS/MS Buffer
9. LC/MS/MS (liquid-chromatography/ tandemmass spectrometry) analysis:

LC- Column: Isocratic using a Hamilton PRP-X200 column. Mass Spectroscopy: MRM mode with parent/daughter ion pairs of 114/44 for creatinine and 117/47 for d3-creatine. (Ref: Takahashi N, Boysen G, Li F, Li Y, Swenberg JA. Tandem mass spectrometry measurements of creatinine in mouse plasma and urine for determining glomerular filtration rate. *Kidney Int.* 71: 266-271, 2006.)

Reagent Preparation:

Reagent 1: Acetonitrile: As supplied by vendor

Reagent 2: LC/MS/MS Buffer: 85% water and 15% of 8.5mM Ammonium Acetate.

Reagent 3: D3-Creatinine standard: (5 mg/dL)