



# Albumin

Version: 1

Edited by: Gary Cline, John Stack: Yale MMPC Analytical Core

**Summary:** Procedure used to determine the concentration of albumin in blood, plasma, and serum. Albumin is measured as its conjugate with bromocresol green monitored at 600 nm.

## Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Albumin Standard	Prolabs	R85260
Albumin Reagent	Prolabs	R85211
Assayed Control Serum 1	Prolabs	R83082
Assayed Control Serum 2	Prolabs	R83083

**Protocol:** Analysis by automated system Cobas Mira Plus.

- 1) Calibrate Cobas for Albumin analysis by running an albumin standard, assayed control serum 1 and assayed control serum 2.
- 2) Sample handling as performed by the Cobas Mira Plus.
  - a) Pipette 2 $\mu$ L of sample into a cuvette slot.
  - b) Add 250  $\mu$ L of Albumin reagent and mix.
  - c) Mixture is incubated at 37°C and spun for 10 minutes.
  - d) Absorbance is measured at 600 nm.

## Reagent Preparation:

Albumin Standard: As supplied by vendor

Albumin Reagent: As supplied by vendor

Assayed Control Serum 1: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.

Assayed Control Serum 2: Add the appropriate amount of water (6.5mL) to the chemical control bottle. Invert to mix, allowing 15 minutes for the reagent to settle.