



β -Hydroxybutyrate

Version: 1

Edited by: Gary Cline, John Stack: Yale MMPC Analytical Core

Summary: Procedure used to determine the concentration of β -Hydroxybutyrate in blood, serum, and plasma. β -Hydroxybutyrate is measured by the oxidation to acetoacetate, followed by reduction of an indicator dye (monitored at 505 nm) by NADH.

Reagents and Materials:

Reagent/Material	Vendor	Stock Number
β -Hydroxybutyrate Controls	Pointe Scientific Inc.	H7587-CTL
β -Hydroxybutyrate Reagent (test kit includes standard).	Pointe Scientific Inc.	H7587-58

Protocol: Analysis by automated system Cobas Mira Plus.

- 1) Calibrate Cobas for β -Hydroxybutyrate by running a β -Hydroxybutyrate standard and three β -Hydroxybutyrate controls.
- 2) Sample handling as performed by the Cobas Mira Plus.
 - a) Pipette 3 μ L of sample into cuvette.
 - b) Add 105 μ L of β -Hydroxybutyrate Reagent.
 - c) Mixture is incubated at 37°C for 10 minutes.
 - d) Absorbance is measured at 505 nm.

Reagent Preparation:

β -Hydroxybutyrate Controls: As supplied by vendor

β -Hydroxybutyrate Reagent and Standard: As supplied by vendor