



SOP: CASE MMPC

Lipid Analysis Assay

by GC-mass spectrometry

#CA2016 / CA2038 / CA2040

Summary:

A known quantity of tissue / plasma is hydrolyzed and extracted after adding known amounts of internal standards: eg. heptadecanoic acid and cholesterol-d₇. Fatty acids / cholesterol are analyzed as their trimethylsilyl derivatives using gas chromatography-electron impact ionization mass spectrometry (GCMS) (note: this protocol outlines the processing for palmitate and cholesterol; other fatty acids and sterols can be assayed using this preparation, see refs 1,2).

Reagents/Materials:

Reagent/Material	Quantity Required	Vendor
KOH / Ethanol	1mL	stock
Heptadecanoic acid and cholesterol-d ₇	25µL 25µL	Sigma Aldrich
HCl	50µL	stock
Chloroform	300µL	stock
*TMS	60µL	Regis

*bis(trimethylsilyl) trifluoroacetamide+ 1% trimethylchlorosilane
(Regis, Morton Grove, IL) (TMS)

Protocol:

- For total bound lipids/cholesterol follow steps 1-12
 - For free lipids/cholesterol weigh tissue as in step 1 (use 1ml HCL 6N in place of KOH ethanol solution- omit step 4, *do not heat*), then proceed to step 6-12
1. Pipette 1ml KOH ethanol solution (1N KOH in 70% EtOH) for every 100mg of tissue or 100 µl of plasma use glass screw top tubes (may use less tissue/plasma)

2. Internals standards (IS): add 25 μ l of 1mg/ml heptadecanoic acid (C17:0) and cholesterol-d₇ for every 100 μ l of tissue (note: adjust added amount of internal standards by testing a representative sample)
3. Homogenize on ice with polytron homogenizer (tissue only)
4. Cover and heat for 3 hours at 85°C
5. Pipette 50-100 μ l of solution into an Eppendorf tube
6. Add 50 μ l of 6N HCl
7. Add 300 μ l of Chloroform
8. Vortex and centrifuge for 2 minutes
9. Take 200 μ l of chloroform phase (bottom layer) and dry in GC vial at 75°C
10. React with 60 μ l of TMS, cover, heat at 75°C for 20 minutes
11. Transfer to GC insert and cap

GC-MS Analysis: Lipid TMS derivatives are analyzed using an Agilent 5973N-MSD equipped with an Agilent 6890 GC system, and a DB-17MS capillary column (30 m x 0.25 mm x 0.25 μ m). The mass spectrometer is operated in the electron impact mode (EI; 70 eV).

12. Selective ion monitoring of mass-to-charge ratios (m/z)
 - a. Palmitate: M0=313-M⁺; IS, C:17=327
 - b. Cholesterol: M0=368-M⁺; IS, Cholesterol-d₇=375
 - c. For other lipids: see references

References:

1. Triglyceride synthesis in epididymal adipose tissue: contribution of glucose and non-glucose carbon sources. Bederman IR, Foy S, Chandramouli V, Alexander JC, Previs SF. J Biol Chem. 2009, 284(10):6101-8.
2. Influence of diet on the modeling of adipose tissue triglycerides during growth. Brunengraber DZ, McCabe BJ, Kasumov T, Alexander JC, Chandramouli V, Previs SF. Am J Physiol Endocrinol Metab. 2003, 285(4):E917-25.