



Triglyceride Protocol

Version: 1
 Replaced by version
 Edited by: Peter Havel - UC Davis Metabolism and Endocrinology Core

- [Summary](#)
- [Reagents and Materials](#)
- [Protocol](#)
- [Reagent Preparation](#)

Summary:

Triglycerides are enzymatically hydrolyzed by lipase to free fatty acids and glycerol. The glycerol is phosphorylated by adenosine triphosphate (ATP) with glycerol kinase (GK) to produce glycerol-3-phosphate and adenosine diphosphate. Glycerol-3-phosphate is oxidized by dihydroxyacetone phosphate (DAP) by glycerolphosphate oxidase producing hydrogen peroxide (H₂O₂). In a Trinder5 type color reaction catalyzed by peroxidase, the H₂O₂ reacts with 4-aminoantipyrine (4-AAP) and 3,5-dichloro-2-hydroxybenzene sulfonate (DHBS) to produce a red colored dye. The absorbance of this dye is proportional to the concentration of triglycerides present in the sample.

Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Calibrator	Fisher Diagnostics	TR43002
Reagents	Fisher Diagnostics	TR22203
PBS		
Microplate		
Platereader		

Protocol:

1. Reconstitute powdered reagent with only 25 ml of distilled water to make a 2X solution.
2. Add 3 µl of calibrator and sample to each well.
3. Add 150 µl of PBS to each well. Read at 540 nm.

IMPORTANT: Make sure not to add any bubbles to the wells when dispensing reagents, this will interfere with reading in the platereader.

4. Add 150 μ l of 2X reagent to each well. Incubate at 37°C for 5 minutes. Read at 540 nm.
5. Subtract blank readings from final readings. The assay will be linear so the unknown samples can be calculated as $(\text{Sample Absorbance} \div \text{Calibrator Absorbance}) \times \text{Calibrator Concentration}$.

Reagent Preparation:

PBS – ready to use

Reagent – reconstitute with distilled water to make a 2X solution